Gebrauchsanweisung
Instruction for Use
Mode d'emploi

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Please read all instructions carefully before using the product. Please comply with all warnings and precautions included in these instructions. Non-compliance may result in serious complications.

1.0 Product description

The conventional treatment of thoracic aortic aneurysm (TAA) and TA dissections is a highly invasive surgical procedure. As is always the case with major surgery, the level of risk is high.

The endoluminal stentgraft prostheses from JOTEC[®] GmbH represent an alternative to the conventional surgical methods and therapeutic practices otherwise applied in treatment of vascular aneurysms, dissections and aortic ruptures of the thoracic aorta.

The JOTEC E-vita[®] thoracic stentgraft system was specifically developed and designed for endoluminal treatment of lesions, in particular aneurysms and dissections in the thoracic aorta descendens.

The E-vita[®] thoracic system has two main components:

- the self-expandable nitinol stent with polyester mantle and platinum-iridium radiopaque markers and
- the delivery system carrying the E-vita[®] thoracic Stentgraft.

The delivery system with the E-vita[®] thoracic Stentgraft is introduced into the body using a 20F, 22F or 24F E-asy sheath or directly via an 0.035" (0.89 mm) rigid, coated guidewire, then advanced up to the damaged section of the thoracic aorta. Retracting the outer shell of the delivery system and separate release of the proximal stent wire activates the self-expandable mechanism of the E-vita[®] thoracic Stentgraft.

The E-vita[®] thoracic Stentgraft covers the damaged vessel segment surface so that the blood flow passes through the lumen of the E-vita[®] thoracic Stentgraft and is shielded from contact with the lesion. If necessary, the E-vita[®] thoracic Stentgraft can be modelled step by step onto the vascular wall with the help of a separate balloon catheter.

1.1 E-vita[®] thoracic TAA Stentgraft

The TAA stentgraft (**Fig. 1 and 2**) comprises individual stent wires (NiTi alloy), permanently sewn into a textile tube (PET) using surgical suture material (PET, woven). E-vita[®] thoracic Stentgrafts can be supplied either in tubular or conical forms, which are realized by varying the diameter and form of the stent wires, providing for optimized adaptation to the damaged vascular segment.

E-vita[®] thoracic Stentgrafts, which have to be implanted far into the aortic arch, feature free wires (Straight Open Design) to hold open the supraaortal vessels, in particular the subclavian artery.

A free wire at the distal end part of the prosthesis holds efferent vessels open. In particularly long lesions, two or more E-vita[®] thoracic Stentgrafts can be implanted so as to overlap. In such cases, the diameters of two adjacent stentgrafts should differ by approximately 20%, at any rate not by more than 6 mm, and the overlap should be at least 30 mm long.

Small radiopaque marker rings made of a Pt-Ir alloy have been sewn into the polyester textile to improve X-ray visibility for positioning and subsequent checking of E-vita[®] thoracic stentgrafts. In the Straight Open – Straight Cut (SO) configuration, 5 radiopaque marker rings are attached and in the Twin Stent – Free Spring / Straight Cut (ST) configuration 3 / 4 radiopaque marker rings are attached.


Figure 1: E-vita[®] thoracic Stentgraft, e.g. article 70SOxxySzz-00


Figure 2: E-vita[®] thoracic Stentgraft, e.g. a) article 70STxxyyFzz-00 b) article 70STxxyySzz-00

1.2 Delivery system

The delivery system (**Fig. 3**) consists of four nesting catheters:

- an internal guide wire catheter
- an intermediate release catheter for proximal stentgraft release
- a pusher (catheter)
- and the outer delivery sheath with proximal textile mantle.

Important: The switch on the release grip must remain in the P position as designated until the stentgraft is released.

The lumen of the guide wire catheter is for the guide wire. The proximal end of the guide wire catheter is fitted with a soft, radiopaque polyurethane (PU) tip. The E-vita[®] thoracic Stentgraft is loaded at the proximal end of the delivery system. The pusher (catheter) is in contact with the distal end of the E-vita[®] thoracic Stentgraft. Shifting the switch on the release handle to position D switches the system to operating mode and the stentgraft can be released. This is realized by means of controlled retraction of the delivery sheath against the pusher. Fixation of the first stent wire can be released at any time during the process by pushing the actuation button (red) on the pusher handle.

Following complete release of the stentgraft, the delivery system is removed from the vascular system.

JOTEC Instruction for use E-vita[®] thoracic


Figure 3: E-vita[®] thoracic Stentgraft System

JOTEC Instruction for use E-vita[®] thoracic

2.0 Indications

The E-vita[®] thoracic Stentgraft from JOTEC[®] GmbH has been specially developed for treatment of aneurysms, dissections and specific lesions of the thoracic aorta descendens. The main indications for use of a TAA stentgraft from JOTEC[®] GmbH would be, in particular, acutely life-threatening patient conditions due to

- aneurysm rupture
- traumatic aortic rupture
- aortic dissection
- a plaque rupture
- a penetrating aortic ulcer (PAU)
- an intramural haematoma (IMH)
- an aortobronchial fistula or aortooesophageal fistula.

The stentgraft was designed for patients for whom conventional surgery is either indicated or contraindicated due to pre-existent risk factors.

3.0 Contraindications

The E-vita[®] thoracic Stentgraft system is contraindicated in the following patients:

- Patients with highly tortuous vessels; contact JOTEC if in doubt.
- Patients with blood vessels or aneurysms not suitable for treatment with the E-vita[®] thoracic Stentgraft system
- Patients in whom the aneurysm or aortic disease to be repaired by the stentgraft includes vitally important vascular exits (carotid arteries, visceral and renal arteries).
- Patients with stenoses that cannot be crossed by the delivery system.
- Patients with a heavy curve on the setting off point of the distal end of the E-vita[®] thoracic Stentgraft.
- Patients in whom the arterial approach site is too narrow for implantation of the E-vita[®] thoracic Stentgraft System.
- Patients in whom the necessary proximal positioning in the aorta is impossible due to the work length of the delivery system.
- Patients in whom the required implantation materials (see below) cannot be used.

The therapeutic indication must be very carefully considered in the following cases:

- Patients with systemic or local infections in whom a bacterial infection of the stentgraft would appear to be a possibility
- Patients in whom administration of radiological contrast agents, anticoagulants or heparin is contraindicated
- Patients in whom the E-vita[®] thoracic Stentgraft System would raise the risk of a paraplegia or loss of sensory function.

Important: The final therapeutic decision must be made by the physician in charge of treatment.

4.0 Sizing of the E-vita[®] thoracic Stentgraft System

Before beginning with the procedure, size and properties of the diseased aortic section and access vessels, such as vessel diameter and status (e.g. calcifications and stenoses, intraluminal thrombi) must be determined with the aid of suitable imaging equipment. 3D spiral CT angiography with contrast agent is recommended. The vessel diameters must be measured in the healthy aortic segment. The stentgraft diameter suitable for the native vessel is determined by “oversizing”.

JOTEC Instruction for use E-vita[®] thoracic

The guideline figure is 20%, but not more than 6 mm, for example vessel diameter 24 mm + 20% = 28.8 mm, corresponding to a stentgraft diameter of 28 mm. In dissections oversizing should not be more than 10%. When determining the length, both the proximal and distal anchoring length must be included in the calculation. At least one covered spring must be placed in the anchoring zone. The distal end of the E-vita[®] thoracic Stentgraft must not be positioned directly in a pronounced kink. This can be avoided by means of careful determination of the length of the E-vita[®] thoracic Stentgraft.

Before beginning with the procedure it must be determined whether the working length of the delivery system is sufficient to reach the required proximal position in the aorta.

Important in cases in which several E-vita[®] thoracic Stentgraft Systems are used:

- The overlap between two or more stentgraft systems should be at least 30 mm long. A longer overlap is recommended in angled or curved segments.
- Further, the stentgrafts must have a diameter difference of 20%, but not exceeding 6 mm
- The stentgraft with the larger diameter must always be placed in the stentgraft with the smaller diameter.

Important: The physician in charge of treatment is responsible for selection of the correct sizes.

To accommodate special anatomical situations, the proximal and distal stentgraft connections can also be custom-made to meet customer requirements. Please contact our Customer Service or a JOTEC representative if you have any questions.

5.0 Preparations for the thoracic stentgraft procedure

5.1 Required equipment for implantation of an E-vita[®] thoracic Stentgraft:

- Radiographic image converter with digital angiography function (mobile C-Arm device or stationary angiography system)
- Fluoroscopic visualization and facilities for recording and callup of all images
- Surgical theatre in case open surgery should become necessary.

Important: The implantation procedure should only be carried out by a physician trained in endovascular techniques (cardiac and thoracic surgeon, vascular surgeon, interventional radiologist or cardiologist). Staff and facilities should be available for simple or complex vascular access or a transition to open surgery.

5.2 Materials required for the implantation procedure

- Various different guide wires, length at least 260 cm
- A rigid, coated guide wire measuring 0.035" (0.89 mm) in diameter and 300 cm in length, preferably the JOTEC E[®]-wire
- Various different angiography and angioplasty catheters
- A suitable sized JOTEC E-asy delivery sheath

Note: Using the JOTEC E-asy delivery sheath reduces the working length of the E-vita[®] thoracic Stentgraft System by 90 mm!

- Sterile delivery sheaths: 7, 8 or 9 French (3 French \cong 1 mm).
- Contrast agent
- Puncture needle
- Sterile syringes
- Heparin and heparinized physiological saline solution
- Balloon catheters


6.0 Heparinization and inhibition of thrombocyte aggregation

To reduce the risk of a thromboembolism, the patient's blood should be anticoagulated during the procedure. As in vascular surgery, between 5,000 and 10,000 IU of heparin, depending on body weight, should be applied intraoperatively. In cases of heparin intolerance, preparations such as hirudin, for example, can be used. Postoperative antithrombocyte therapy may be indicated as judged by the physician in charge of treatment.

Caution: *Insufficient heparinization may result in thrombosis and embolism.*

7.0 Preparation for implantation

- The access vessel, a. iliaca or a. femoralis, is exposed using standard surgical techniques.
- A second point of access can be used for diagnostic purposes to lay an angiography catheter. This access point can be located as per the judgement of the physician percutaneous, femoral or brachial, etc.
- Under fluoroscopy, a fitting E-asy sheath for the delivery system can be inserted up to immediately below the renal arteries using an 0.035" (0.89 mm) guidewire.
- Insert an angiography catheter through the second access point to localize the implantation site. The angiography catheter can stay put during the implantation procedure for orientation purposes.
- After removing the two protective wires from the delivery system, the lumen of the guide wire catheter is flushed with heparinized saline solution.
- Make sure the flushing tube is covering the textile tube in which the E-vita[®] thoracic Stentgraft is loaded. The delivery system is held with the proximal end up and is flushed with heparinized saline solution through the inserted flushing valves at the pusher handle.

Important: *The E-vita[®] thoracic must not be bent or kinked without a guide wire; make sure the catheters have not been damaged. In case of doubt, the system must be discarded and replaced by an intact system.*


7.1 Delivery of the E-vita[®] thoracic Stentgraft System

- The E-vita[®] thoracic Stentgraft System is delivered into the vessel by means of a rigid, coated 0.035" (0.89 mm) guide wire. Use a matching E-asy sheath in case of calcified and/or stenosed vessels.
- The E-vita[®] thoracic Stentgraft must be completely deaerated just before it is introduced into the artery or delivery sheath. If necessary, introduce a further bolus of flushing liquid (heparinized saline solution) into the stentgraft through the inserted flushing valve on the pusher handle. The flushing tube must cover the loaded stentgraft. Then introduce the E-vita[®] thoracic Stentgraft System into the exposed artery or in the indwelling sheath, whereby the flushing tube is pushed back of its own accord.
- After the textile sheath is completely inserted into the access vessel or the sheath remove this flushing tube.
- The E-vita[®] thoracic Stentgraft System is slowly advanced to the implantation site under fluoroscopy.
- When the pusher handle needs to be rotated, use fluoroscopy to make sure, tip and stent follow the rotation. Do not twist the catheter system as this could lead to catheter damage!
- Prior to release of the E-vita[®] thoracic Stentgraft, the exact position and situation of the stentgraft must be checked by means of angiography. Then fix the image converter and mark the orientation points for release (roadmap, mark on screen or external marking, e.g. with needles). The E-vita[®] thoracic Stentgraft features two radiopaque marking rings at the proximal end (**Fig. 1**). These markers form the beginning of the circular graft mantle. Further markers are located in the middle and two others at the distal end of the stentgraft.

JOTEC Instruction for use E-vita[®] thoracic

- If the E-vita[®] thoracic Stentgraft is to be located on an efferent vessel, e.g. the a. subclavia, the first radiopaque ring must be positioned **distally** in relation to the efferent vessel.
- As soon as the E-vita[®] thoracic Stentgraft System has been precisely positioned, it may be appropriate to reduce the patient's mean arterial blood pressure to about 60 mmHg. This measure is taken to prevent unintentional shifting of the E-vita[®] thoracic Stentgraft upon release.
- If still necessary, the transport retainer should now be removed (white Velcro band).
- Push the switch on the release handle to position D.
- The orange lever on the handle is raised. The system is ready for release.

Release handle
opened (Position D)


- By raising and lowering the orange lever on the release handle while retaining the pusher handle with the other hand the stentgraft is released from proximal to distal.

Important: The position of tip and stentgraft is controlled with the pusher handle.

- The textile sheath with the marking ring is retracted until the first three stent wires are released (Fig. 4). This is checked under fluoroscopy. Thereby the first wire is still held proximally at the tip of the guide wire catheter.
- The position of the E-vita[®] thoracic Stentgraft in relation to the implantation site should be checked once again using angiography or fluoroscopy. If the E-vita[®] thoracic Stentgraft has partially opened too far below or above the intended implantation site, the entire delivery system can be carefully pushed in distal or proximal direction until the first X-ray marker shows the desired position.

Important: This manoeuvre cannot be performed if the proximal stent wire has been released unintentionally!

- The proximal stent wire can now be released from the delivery system by removing the white safety retainer from the pusher handle, then pushing the red actuation button. Check under fluoroscopy to make sure the proximal wire is completely released.

Important: Replace the white safety retainer at the intended location on the pusher handle once release is realized.

- The E-vita[®] thoracic Stentgraft is completely released in distal direction until the release handle contacts the fixated pusher handle.
- Finally, the switch on the release handle is pushed to position N. Now fixate the release handle and pull the delivery system back into the initial positions using the pusher handle. Then remove the closed delivery system from the aorta and sheath, respectively.
- The guidewire must not be pulled back or removed when removing the E-vita thoracic delivery system.
- For treatment of longer aortic segments, additional E-vita[®] thoracic Stentgrafts can be used. Ensure overlapping of at least 30 mm and oversizing of 20%, although not more than 6 mm. The stentgraft with the larger diameter must be positioned / opened in the stentgraft with the smaller diameter.

JOTEC Instruction for use E-vita[®] thoracic

- Having excluded an endoleakage and documented the complete exclusion of the aortic lesion in the final angiography, the delivery sheath, if used, is removed along with the stiff guide wire.
- Direct intraoperative control of therapeutic success, i.e. exclusion of leakage and patency of the stentgraft, is best realized by means of digital subtraction angiography (DSA).
- The arterial access point is closed using standard surgical techniques.

8.0 Warnings and precautions

- Do not use the E-vita[®] thoracic Stentgraft System with patients in whom the necessary preoperative and postoperative imaging examinations cannot be carried out. All patients should be carefully monitored and must undergo regular checkups to determine changes in their pathological status and the functionality and integrity of the stentgraft.
- Do not use the E-vita thoracic Stentgraft if the work length of the delivery system is insufficient.
- Do not use the E-vita thoracic Stentgraft System if the aorta is heavily curved on the setting off point of the distal end.
- This product must not be resterilized. The product must not be used if the packaging is damaged or has been opened outside of the sterile area.
- Do not use the E-vita thoracic Stentgraft System if the proximal end of the textile tube has slipped off the tip.
- Do not use the E-vita thoracic Stentgraft System if the proximal or distal end is bent or buckled (possible transport damage).
- Insufficient heparinization may result in a thrombosis or embolism.
- Patient sensitivity or tendency to show an allergic reaction to heparin and/or contrast agent must be kept in mind.
- The E-vita[®] thoracic Delivery System must be disposed of in accordance with valid hospital, administrative and other official regulations.
- The E-vita[®] thoracic Stentgraft in the delivery system must be completely deaerated (ventilated) shortly before it is inserted into the arterial vascular system. If necessary, inject a bolus of flushing liquid (heparinized saline solution) into the stentgraft through the inlet valve on the pusher handle with the flushing tube pulled over it. The textile mantle of the delivery system ensures complete deaeration, even in a horizontal position.
- The E-vita[®] thoracic Stentgraft must not be twisted prior to release.
- The switch on the release handle must remain in position P until the stentgraft is released.
- Rotational manoeuvres to position the proximal E-vita[®] thoracic Stentgraft must not be carried out in the aortic arch or in other highly curved aortic segments.
- Only the mantled segment of the E-vita[®] thoracic Stentgraft may be modelled to the aortic wall with the balloon. Modelling of a freestanding, non-mantled wire may result in perforation of the aortic wall.
- When retracting the guide wire, make sure the distal part of the tip does not become entangled in the proximal stent wire. Otherwise, the E-vita[®] thoracic Stentgraft could be dislocated or collapsed. Please observe the position of the stentgraft carefully during this manoeuvre. Review the situation if resistance is felt during the retraction, always keep the position of the guidewire under observation.
- **Leaks** are the single most frequently observed complication in endovascular stentgraft implantations. Following implantation of the E-vita[®] thoracic Stentgraft it is important to check for leakage at the stentgraft in the aneurysm or in the false dissection lumen. Minimal, temporary leaks through the prosthesis textile may be due to the intraoperative heparinization. Checks for leaks must be included in follow-up examinations.
- **Microembolisms** – The risk of microembolism rises with the amount of manipulation and the duration of the procedure, so that excessively manipulative and long endovascular procedures should be avoided.
- **Vascular thrombosis** – Folds and wrinkles in the stentgraft may cause thrombosis formation. Subsequent dilatation of the E-vita[®] thoracic Stentgraft with a separate balloon should therefore be practised to reduce folds and wrinkles to a minimum. Always make sure the balloon does not expand or blow up too much. Excessive stenosing of the E-vita[®] thoracic Stentgrafts may also lead to formation of thromboses. Following the stentgraft placement an intraoperative angiogram or postoperative chest x-ray must ensure that the implanted E-vita[®] thoracic Stentgraft is neither buckled nor stenosed. Highly stenosed (narrowed) segments require immediate treatment.
- **Migration** – If the E-vita[®] thoracic Stentgraft is too small, distal stent migration is a possibility. When sizing the E-vita[®] thoracic Stentgraft 20% oversizing (but no more than 6 mm), in case of dissections no more than 10% is recommended. Other causes of migration include insufficient fixation within the total available aortic neck length and

JOTEC Instruction for use E-vita[®] thoracic

expansion of the proximal stentgraft in a vessel that is filled with thrombi or characterized by a sharp bend. The physician in charge of treatment is responsible for selecting the proper dimensions.

- **Overdilatation** – Overstretching and damage may develop if the E-vita[®] thoracic Stentgraft size used is too large in relation to the vessel diameter; such damage can also result from excessive expansion/modelling pressure applied by means of a subsequent balloon manoeuvre.
- **Stroke** – Embolism fragments or blockage of the brachiocephalic vessels may cause a stroke.
- **Intestinal ischaemia** – Complete or partial occlusion of the tr. coeliacus and a. mesenterica superior by an incorrectly positioned E-vita[®] thoracic Stentgraft or stent-caused embolisms may lead to ischaemia of the organs of the upper abdomen and intestine.
- **Renal insufficiency** – Complete or partial occlusion of the a. renalis by an incorrectly positioned E-vita[®] thoracic Stentgraft or stent-caused embolisms may lead to severe renal complications. Excessive amounts of contrast agent may also cause serious renal complications.

9.0 Possible complications and side effects

The literature reports 30-day mortality rates of between 7 and 10% following stentgraft implantations. The 30-day mortality rates for the alternative procedure, open surgery, are between 12 and 31%. Complications may occur at any time during or after the stentgraft implantation. Physicians intending to perform the procedure must be familiar with the possible complications. Late complications may also occur and the course should be monitored in regular aftercare checkups.

Clinical complications:

- Potential vascular complications include: thrombosis, thromboembolism, occlusion (arterial and venous), dissections, perforation, occlusion of collateral vessels, vascular ischaemia, tissue necrosis, amputation
- Vascular rupture / aneurysm
- Potential neurological complications include: paralysis, paraparesis, cerebrovascular events, transitory ischaemic episodes and stroke.
- Intestinal ischaemia
- Infection and fever
- Gastrointestinal complications
- Haemorrhaging
- Pulmonary complications
- Wound healing complications
- Oedema
- Cardiac failure / infarction
- Transition to open surgery
- Death

Problems involving the E-vita[®] thoracic Stentgraft System:

- Endoleaks
- Migration
- Jamming, entanglement of various delivery system components (release blockade)
- Disassembly of release handle from the outer introducer sheath in case of release blockade
- Twisting or buckling of the stentgraft
- Difficulties with delivery or removal of the delivery system
- Expansion difficulties (e.g. stentgraft gets stuck in delivery system)
- Stentgraft does not open up
- Inexact positioning (e.g. false orientation of the stentgraft)
- Tearing of the graft material
- Breakage of stentgraft wires

JOTEC Instruction for use E-vita® thoracic

10.0 MRI Information

Safety information concerning magnetic resonance imaging (MRI) refer to use of shielded MRI systems with static magnetic fields of 1.5 Tesla or less, gradient magnetic fields of 20 Tesla/second or less and a maximum mean total body specific absorption rate (SAR) of 1.2 W/kg during a 30-minute imaging session.

The stentgraft is MRI-safe, that is the stentgraft in the body of a patient undergoing an MRI procedure does not represent an additional risk, although it may influence the MRI image quality depending on the pulse sequence and the area to be evaluated.

11.0 Storage and packaging

The E-vita® thoracic Stentgraft System is a single-use (disposable) product sterilized by gamma radiation which is supplied in a sterile package and remains sterile as long as the package has not been opened or damaged.

The stentgraft system must be stored in a dark, dry, cool location at temperatures of between 10°C and 38°C. Avoid direct exposure to sunlight, excessive heat can have considerable negative effects on the functionality of the product. Do not expose the stentgraft system to ionizing radiation or ultraviolet light.

Guarantee and warranty

JOTEC® GmbH guarantees that this product complies with the legal requirements, in particular the EU Directive on Medical Products 93/42/EEC. JOTEC® GmbH will therefore, if the product demonstrably has material or manufacturing defects prior to its expiry date that have resulted during the manufacturing or packaging procedures:

- a) replace the product with a product of the same type and design or
- b) with a product that is of equal value as to functionality.

All of the following conditions must be met before the guarantee applies:


- a) The product must not have been used after its expiry date.
- b) The defect did not occur for the first time after the product expiry date.
- c) For products approved for single use only or intended by JOTEC® GmbH for single use only, product use revealing the defect must have occurred during first use.
- d) The defect must be reported to JOTEC® GmbH or its representative without delay.
- f) The defective product must have been sent to JOTEC® GmbH or its representative immediately following reporting of the defect.

The above guarantee does not affect any legal regulations concerning warranty liability by the seller or provisions of product liability law and applies in addition to such regulations.

Otherwise, the guarantee applies only in the scope described above. Further guarantee claims are excluded. The guarantee and the legal claims mentioned above apply in place of and to the exclusion of any other provisions by JOTEC® GmbH, whether made orally, in writing, expressly or tacitly or in any other manner. None of our representatives or employees has the right to establish any other provisions or make agreements to contrary effect.

Due to many factors beyond the influence of JOTEC® GmbH such as shipping, storage, handling by users, misuse and non-compliance with instructions for use or warnings, JOTEC® GmbH cannot be made liable for damages of any kind caused directly or indirectly by the product or its use if the damage is not demonstrably caused by product defects. JOTEC® GmbH cannot assume liability for false sizing of the stentgraft, improper use of the product or false positioning of the E-vita® thoracic Stentgraft. JOTEC® GmbH also cannot be made liable in particular for damages caused by resterilization or reuse of products approved or intended by JOTEC® GmbH for single use only.

Appendix:

1. Prepared E-vita[®] thoracic System
2. Opening of the release mechanism by pushing the switch on the release handle from P to D.
3. Start of release by means of steady raising and lowering the opened orange lever. The pusher handle is held in place.
4. Release of the proximal fixation after the white retainer has been removed from the pusher handle and the red actuation button has been pushed. Afterwards, the white retainer is returned to its former position on the pusher handle.
5. Complete release of stentgraft.

Figure 4: Release procedure for E-vita[®] thoracic Stentgraft, steps 1 to 5

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JOTEC GmbH